

# Protocol 001 - Preparing RT Collagen Coated Flasks for Chordoma Cell Line Culture

---

## Introduction

This protocol will cover coating T flasks with Rat Tail Collagen I (RT Col I) for general chordoma cell line culture.

## Materials

- > [Sigma Rat Tail Collagen \(Sigma, Cat# C3867-1VL\)](#) **\*\*\*CHECK LOT NUMBER FOR CONCENTRATION\*\*\***
- > Sterile ddH<sub>2</sub>O
- > T flasks
- > 1X PBS, pH 7.4 (Thermo Fisher, Gibco, Cat# 10010023)

## Procedure

### Cell Culture Flask Preparation - Thawing Procedure

1. **\*\*\*CONCENTRATION IS LOT DEPENDENT. USE "SOURCE" ON VIAL AS LOT NUMBER. CONCENTRATION TO BE WRITTEN ON 15mL CONICAL ALIQUOTS\*\*\*** Dilute 2.9mg/mL of Rat Tail Collagen I at a **1:58 ratio** in ddH<sub>2</sub>O for a final concentration of **50 ug/mL**  
  
**\*\*\*RT Collagen may have solidified in 4C. If this is the case just warm the collagen to room temperature, in the dark, and wait for it to melt\*\*\***
2. For **5x T25** flasks dilute **224 uL** of RT Col I into **12.78 mL** of sterile ddH<sub>2</sub>O.  
  
For **5x T75** flasks, multiply this by 3x. (3\*224uL of RT Col I and 3\*12.78mL of Sterile ddH<sub>2</sub>O = **672 uL** of RT Col I and **38.34 mL** of sterile ddH<sub>2</sub>O)
3. Add **2.5 mL** of RT Col I dilution to **T25** flasks and **7.5 mL** to **T75** flasks
4. Incubate at **37°C** for **1 hour**
5. Aspirate Rat Tail Collagen I Dilution and thoroughly wash with equal volume of 1X PBS, pH 7.4. Repeat 2x
6. Immediately use flasks or store at **4°C** for up to **1 month**