

Protocol 013 - Freezing Chordoma Cell Lines (general)

Introduction

This is a general protocol for freezing chordoma cells lines

Materials

- > Cell culture media
- > FBS
- > DMSO (filtered sterile)
- > 2mL Cryovial
- > 15mL/50mL conical
- > 0.2um nitrocellulose filter
- > Appropriate size syringe for filter
- > Mr. Frosty with fresh isopropanol

Procedure

Freezing Procedure

1. **CRITICAL** Cells should be frozen when they are ~80% - 90% confluent in a vessel that was dedicated for freezing. The cells should also be tested for mycoplasma to circumvent freezing contaminated lines
2. Make 20% FBS and 10% DMSO freezing media
 - a. One T75 flask = three cryovials = 3 mL of freezing media (1 mL / cryovial)
 - b. Since we are filtering and volume will be lost, we will need to make a couple mL's extra
 - c. So, one 5 mL of freezing media = one T75 flask = 4 mL of 10% FBS media + 500 uL of FBS + 500 uL of DMSO (sterile filtered)
 - CRITICAL** d. Filter the freezing media sterile with a 0.2um nitrocellulose filter *****use nitrocellulose filter for DMSO compatibility*****
3. Label the cryovials with:
 - cell line name
 - passage number
 - % confluence when thawed in a T75 flask (so 30% if three cryovials were made from a single 90% confluent T75 flask)
 - date
4. Aspirate media from 80% - 90% confluent flask

5. Wash cells with appropriate amount of PBS (8 mL / 75 cm²) one time and aspirate PBS fully from flask
6. Add appropriate amount of trypsin (3 mL / 75 cm²) and place flask in the incubator at 37°C and 5% CO₂ until cells have fully detached
7. After cells have fully detached, resuspend them in an appropriate volume of media (8 mL / 75 cm²)
8. Spin the cells at 250 rcf for 2.5 min
9. Carefully aspirate the supernatant without disrupting the cell pellet
10. Resuspend the cell pellet in 3 mL of freezing media and aliquot 1 mL of the suspension to each cryovial, so three cryovials in total
11. Place the cryovials in the Mr. Frosty and store overnight in the -80°C
12. The next day remove the cryovials with frozen cell suspension to the LiN2 for long term storage