

CH22 Growth Conditions

Growth Conditions

Media:


RPMI + 10% FBS + Pen/Strep

→ Passage when ~80-90% confluent (1:2, 1:3)


→ Change media every 2-3 days

Phase Contrast Image Review

One vial of cells arrived. Were thawed on 1/6/2015. Thaw produced considerable 24 hour cell death but culture was visibly clear of contamination and recovered well after media change.



24 Hours Post Thaw, Pre-media change (10/08/2013)



24 Hours Post Thaw, Post-media change (10/08/2013)